Pyonpyon3 Behavioral Program:

Kentaroh Takagaki, M.D. Ph.D Updated 2015.10.26

Inputs

Outputs

Sound

Sound AO: Audio output (if output from DAQ card)

Sound card: Sound card number (if output from computer sound card—up to ?10ms jitter)

Shock Intensity

Shock intensity AO:

Code/LED DO (and Tone gate PFI)

Ī	Trial Type Number								Current State				Null	Sw	Tone	Stim
ĺ	15	14	13	12	11	10	9	8	7	6	5	4	3	2	1	0

Bit 0: Stimulation gate (LED output)

Bit 1: Tone gate (LED output)

[Tone gate PFI: PFI channel same as bit 1 above, used to trigger sound analog output]

Bit 2: Sw1/animal left (LED output)

Bit 3: Reward


Bit 4-7: Current trial state

Bit 8-15: Current stimulus trial type

Switches DI

Webcam

Flow Chart



Preparation code